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鼻咽頭抗原定性 : ESPLINE SARS-CoV-2 rapid antigen test (Fujirebio)

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TABLE 2. Test results and performance characteristics of the Abbott BinaxNOW COVID-19 Ag Card Point of Care Diagnostic Test (BinaxNOW antigen test) compared with real-time reverse transcription–polymerase chain reaction (RT-PCR) for testing received among asymptomatic and symptomatic persons at two community-based testing sites — Pima County, Arizona, November 2020

Results and Performance	Real-time RT-PCR, no. of tests		
	Positive	Negative	Total
<b>BinaxNOW antigen test result</b>			
<b>All participants (N = 3,419)</b>			
Positive	157	4	161
Negative	142	3,116	3,258
<b>Total</b>	<b>299</b>	<b>3,120</b>	<b>3,419</b>
<b>Symptomatic (≥1 symptom) (n = 827)</b>			
Positive	113	0	113
Negative	63	651	714
<b>Total</b>	<b>176</b>	<b>651</b>	<b>827</b>
<b>Asymptomatic (n = 2,592)</b>			
Positive	44	4	48
Negative	79	2,465	2,544
<b>Total</b>	<b>123</b>	<b>2,469</b>	<b>2,592</b>
<b>BinaxNOW antigen test performance, % (95% CI)</b>			
<b>All participants (N = 3,149)</b>			
Sensitivity	52.5 (46.7–58.3)		
Specificity	99.9 (99.7–100.0)		
PPV	97.5 (93.8–99.3)		
NPV	95.6 (94.9–96.3)		
<b>Symptomatic (n = 827)</b>			
Sensitivity	64.2 (56.7–71.3)		
Specificity	100.0 (99.4–100.0)		
PPV	100.0 (96.8–100.0)		
NPV	91.2 (88.8–93.1)		
<b>Asymptomatic (n = 2,592)</b>			
Sensitivity	35.8 (27.3–44.9)		
Specificity	99.8 (99.6–100.0)		
<b>PPV</b>	<b>91.7 (80–7.7)</b>		
NPV	96.9 (96.1–97.5)		


**Abbreviations:** CI = confidence interval; COVID-19 = coronavirus disease 2019; NPV = negative predictive value; PPV = positive predictive value.

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FIGURE. Abbott BinaxNOW COVID-19 Ag Card Point of Care Diagnostic Test (antigen test) results, N1 cycle threshold (Ct) values,\* and viral culture results† among A) symptomatic (N = 136)<sup>‡</sup> and B) asymptomatic (N = 88)<sup>§</sup> participants receiving positive SARS-CoV-2 real-time reverse transcription–polymerase chain reaction (RT-PCR) test results at two community-based testing sites — Pima County, Arizona, November 2020

\* Only those specimens that were analyzed using the CDC 2019-nCoV Real-Time RT-PCR Diagnostic Panel for detection of SARS-CoV-2 and that were analyzed using viral culture are included in the graph.

† Twenty specimens with Ct values <18 had positive antigen and real-time RT-PCR results but were culture negative. The culture showed evidence of cytopathic effects and had presence of SARS-CoV-2 RNA as detected by real-time RT-PCR in the first passage culture, but viral recovery was not two Ct values lower than the corresponding clinical specimen Ct.

‡ Antigen test results: 88 positive and 48 negative; median Ct values indicated with black line: 22.0 for antigen-positive specimens and 33.9 for antigen-negative specimens.

§ Antigen test results: 37 positive and 51 negative; median Ct values indicated with black line: 22.5 for antigen-positive specimens and 33.9 for antigen-negative specimens.

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
TABLE 2. Sensitivity, specificity, positive predictive value, and negative predictive value of Sofia SARS Antigen Fluorescent Immunoassay compared with real-time reverse transcription-polymerase chain reaction (RT-PCR) among asymptomatic and symptomatic persons — two universities, Wisconsin, September–October 2020

Antigen test result	RT-PCR result, no.					
	Asymptomatic (N = 871)			Symptomatic* (N = 227)		
	Positive	Negative	Total	Positive	Negative	Total
Positive	7	14	21	32	2	34
Negative	10	840	850	8	185	193
Total	17	854	871	40	187	227
Test evaluation, % (95% CI)						
Sensitivity	41.2 (18.4–67.1)			80.0 (64.4–90.9)		
Specificity	98.4 (97.3–99.1)			98.9 (96.2–99.9)		
Positive predictive value	33.3 (14.6–57.0)			94.1 (80.3–99.3)		
Negative predictive value	98.8 (97.8–99.4)			95.9 (92.0–98.2)		

Abbreviation: CI = confidence interval.

\* One or more symptoms reported.

FIGURE. Viral culture results among participants with positive Sofia SARS Antigen Fluorescent Immunoassay or positive SARS-CoV-2 real-time reverse transcription-polymerase chain reaction (RT-PCR) results (n = 69),\* by cycle threshold (Ct) value† and the interval between specimen collection and reported symptom onset or asymptomatic status — university A, Wisconsin, September–October 2020

\* n = 30 antigen- and culture-positive; n = 22 antigen-positive and culture-negative; n = 15 antigen- and culture-negative; n = two antigen-negative and culture-positive.

† Ct values represent cycle thresholds for the N1 target probe during SARS-CoV-2 real-time RT-PCR; Ct values are represented on the y-axis in descending order to indicate that lower Ct values represent higher levels of RNA in the specimen.

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
High frequency testing with low analytic sensitivity versus low frequency testing with high analytic sensitivity.

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## **Potential for False Positive Results with Antigen Tests for Rapid Detection of SARS-CoV-2 - Letter to Clinical Laboratory Staff and Health Care Providers**

The U.S. Food and Drug Administration (FDA) is alerting clinical laboratory staff and health care providers that false positive results can occur with antigen tests, including when users do not follow the instructions for use of antigen tests for the rapid detection of SARS-CoV-2. Generally, antigen tests are indicated for the qualitative detection of SARS-CoV-2 antigens in authorized specimen types collected from individuals who are suspected of COVID-19 by their healthcare provider within a certain number of days of symptom onset. The FDA is aware of reports of false positive results associated with antigen tests used in nursing homes and other settings and continues to monitor and evaluate these reports and other available information about device safety and performance.

The FDA reminds clinical laboratory staff and health care providers about the risk of false positive results with all laboratory tests. Laboratories should expect some false positive results to occur even when very accurate tests are used for screening large populations with a low prevalence of infection. Health care providers and clinical laboratory staff can help ensure accurate reporting of test results by following the authorized instructions for use of a test and key steps in the testing process as recommended by the Centers for Disease Control and Prevention (CDC), including routine follow-up testing (reflex testing) with a molecular assay when appropriate, and by considering the expected occurrence of false positive results when interpreting test results in their patient populations.

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**Table 1. Situations where SARS-CoV-2 Ag-RDTs should not be used, based on currently available information**

<b>Do not use SARS-CoV-2 Ag-RDTs:</b>	<b>Explanation</b>
In individuals without symptoms unless the person is a contact of a confirmed case	Pre-test probability (the likelihood, before testing, that the patient has the disease based on epidemiology, case contact, clinical findings) is low.
Where there are zero or only sporadic cases	Ag-RDTs are not recommended for routine surveillance purposes or case management in this setting. Positive test results would likely be false positives. Molecular testing is preferred.
Appropriate biosafety and infection prevention and control measures (IPC) are lacking	To safeguard health workers, respiratory sample collection for any test from patients with suspected COVID-19 requires that operators wear gloves, gown, mask and face shield or goggles (19, 22, 23).
Management of the patient does not change based on the result of the test	If test-positive and test-negative patients will be treated the same way because of unknown or low PPV and/or NPV, then there is no benefit to testing.
For airport or border screening at points of entry	Prevalence of COVID-19 will be highly variable among travellers, and it is therefore not possible to determine PPV and NPV of test results. Positive and negative tests would require confirmatory testing to increase PPV and NPV for decision making.
In screening prior to blood donation	A positive RDT result would not necessarily correlate with presence of viremia. Asymptomatic blood donors do not meet the definition of a suspect case (24).

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